

ul. Kolejowa 33 E, 18-218 Sokoły, tel. 86 274 13 42, fax 86 274 13 50 PKO BP S.A. 96 1020 1332 0000 1502 0324 0447 NIP: 966-176-77-14, REGON: 200133730 www.termopasty.pl, e-mail: biuro@termopasty.pl

Product Data Sheet



Gel Flux

Flux gel is the condensed rosin type flux of RMA class, 1.1.2/3 A type according to ISO 9454-1, intended for SMT assembly and repairs. Flux gel can be applied through a sieve, pattern or from a syringe. Activators utilized in flux manufacturing process enable to leave its remnants after soldering on the package (as far as the process does not include washing phases).

Application:

- SMT assembly
- SMT repair

Package:

Capacity	Package	Symbol
1,4 ml (syringe)	5	ART.AGT-047
14 cm ³ (syringe)	2 or 16	ART.AGT-088
100 cm ³ (plastic can)	1 or 4	ART.AGT-089

